



Four Ways to Optimizing Bioprocessing Bag Reliability, Functionality, and Cost.

Single-use bioprocessing bags have revolutionized the biopharmaceutical industry, offering cost-effective and flexible solutions for the storage and transport of bioproducts. The development of these bags involves meticulous selection of film structures and components to ensure they meet stringent regulatory and performance standards. This paper discusses the challenges associated with selecting the appropriate film and components, the necessity of design for manufacturability assessments, common mistakes in the design process, and strategies to mitigate these issues.

Challenges in Selecting Film Structures and Components

Film Structure Selection

The film structure is critical as it determines the bag's mechanical strength, chemical compatibility, gas permeability, and sterilization capacity. Key considerations include:

Chemical Compatibility: Ensuring the film does not interact negatively with the bioproduct.

Mechanical Strength: Adequate tensile strength and puncture resistance to withstand handling and processing stresses.

Gas Barrier Properties: Minimizing oxygen and CO₂ permeability to maintain product stability.

Sterilization Compatibility: Films must withstand gamma irradiation or other sterilization methods without compromising integrity.

Component Selection

Components such as ports, connectors, and filters must be compatible with the film and the intended bioprocessing application. Key considerations include:

Material Compatibility: Components must be chemically inert and compatible with the film and bioproduct.

Ease of Integration: Components should seamlessly integrate into the bag design to ensure reliability and functionality.

Regulatory Compliance: All materials must comply with regulatory standards such as USP Class VI and ISO 10993.

Design for Manufacturability Assessment

A design for manufacturability (DFM) assessment is critical to identify potential manufacturing challenges early in the design process. This involves evaluating:

Material Availability: Ensuring materials are readily available and have consistent quality.

Manufacturing Processes: Assessing the feasibility of manufacturing processes, including sealing techniques and assembly procedures.

Quality Control: Establishing robust quality control measures to detect defects and ensure consistency.

Common Mistakes

- **Inadequate Film Testing:** Failing to thoroughly test the film for mechanical and chemical properties.
- **Ignoring Component Compatibility:** Overlooking the compatibility of components with the film and bioproduct.
- **Insufficient Quality Control:** Implementing inadequate quality control measures leading to product variability and defects.
- **Incorrect Tooling Design:** Not designing production tools properly, particularly welding tools, which must create a liquid-tight seal. Port welds are high-risk areas where improper tooling design can lead to leaks.



Vonco leads in designing and developing single-use bioprocessing bags through unparalleled expertise in technical materials and design.

Designing and developing single-use bioprocessing bags requires careful consideration of film structures and components to ensure they meet the necessary performance and regulatory standards. Conducting a comprehensive design for manufacturability assessment and addressing common pitfalls through rigorous testing, compatibility assessments, and robust quality control can significantly enhance the reliability and functionality of these critical bioprocessing tools.

THE CONCLUSION

Designing and developing single-use bioprocessing bags requires careful consideration of film structures and components to ensure they meet the necessary performance and regulatory standards. Conducting a comprehensive design for manufacturability assessment and addressing common pitfalls through rigorous testing, compatibility assessments, and robust quality control can significantly enhance the reliability and functionality of these critical bioprocessing tools.

THE VONCO SOLUTION

Vonco Products stands at the forefront of designing and developing single-use bioprocessing bags, offering unparalleled expertise in technical materials and design for manufacturability. Our comprehensive approach leverages decades of experience to ensure the optimal selection of film structures and components, tailored to meet stringent bioprocessing requirements. Utilizing our robust ISO 13485 quality system, we guarantee consistent and reliable product performance through rigorous quality control measures. By partnering with Vonco Products, customers can expedite their design and development processes, leveraging our technical skills to avoid common pitfalls and costly mistakes. Our dedication to innovation and quality ensures that our bioprocessing solutions not only meet but exceed industry standards, providing our clients with a competitive edge in the biopharmaceutical market.